Improving African Swine Fever Surveillance

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

Ha Noi Office: No. 32,
22 Lane, Hung Thinh,
Yen So, Hoang Mai
District, Ha Noi.

HCM Office: No.5 Ta
Hien Street, Thanh My
Loi Ward, Thu Duc City,
Ho Chi Minh City.

**GOLD
STANDARD
DIAGNOSTICS**

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SUPPLIERS

KIT CHẨN ĐOÁN THÚ Y

Lateral Flow

ELISA

PCR

- Viruses
- Bacterial pathologies
- Farm and companion animals

:CZ VACCINES

THIẾT BỊ

- Lab Equipments

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- Micropipette tips
- Tubes

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African Swine Fever Virus

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AFRICAN SWINE FEVER (ASF) Situation Report 64

Period covered: April 2025

This report provides an update of the ASF situation, according to the information shared with WOAHP.

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The spatial distribution of outbreaks reported since January 2022 in domestic pigs and wildlife is shown in Figure 1.


Figure 1. Map of ASF outbreaks which started during 01 Jan 2022 – 31 Mar 2025 in domestic pigs and wildlife.

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RENEWED!

11.AS2.K42: INGEZIM ASF CROM Ag 2.0

Key product for early detection on the field (Point-of-Care). Especially efficient for new outbreaks

RENEWED!

11.ASF.K3: INGEZIM® PPA COMPAC 2.0

Test validated by the European Reference Lab CISA-INIA and recommended by the OIE Reference Lab VISAVET

11.ASF.K1: INGEZIM ASFV-R

Particularly useful for “non-well-preserved” samples

11.PPA.K2: INGEZIM PPA DAS 2.0

Recommended as a herd virus detection, allows large scale sample testing without special lab equipment

RENEWED!

11.ASF.K41: INGEZIM ASF CROM Ab 2.0

Detection in pig and wild boar serum, plasma and blood samples. Used recommended combined with CROM Ag

NEW!

11.ASF.K.43: INGEZIM ASF Combo CROM Ag/Ab



Detection in pig and wild boar serum, plasma and blood samples.


RENEWED!

11.PPA.K.5TX/Q: INGENE q PPA (ASFV)

Real time kit with improved diagnostic performance, better performance than OIE method for pooled testing

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

Evaluation of LFAs on animals infected with Georgia10 (genotype II) ASFV strain

Table 2: Summary results from animals infected with Georgia10 (genotype II) ASFV strain.

Virus quantity (HAD ₅₀)	Clinical scoring	Assay A	Assay B	Assay C	Assay D	Assay E
Day 0						
Animal 1	3	NEG	NEG	NEG	NEG	NEG
Animal 2	4	NEG	NEG	NEG	NEG	NEG
Animal 3	1	NEG	NEG	NEG	NEG	NEG
Animal 4	2	NEG	NEG	NEG	NEG	NEG
Day 4						
Animal 1	9	POS	POS	NEG	POS	POS
Animal 2	8	NEG	POS	NEG	NEG	NEG
Animal 3	1	NEG	NEG	NEG	NEG	NEG
Animal 4	6	NEG	NEG	NEG	NEG	NEG
Day 6 (PM)						
Animal 1	NS	NS	NS	NS	NS	NS
Animal 2	17	POS	POS	POS	POS	POS
Animal 3	16	POS	POS	POS	NEG	POS
Animal 4	11	POS	POS	POS	NEG	POS
Day 8 (PM)						
Animal 1	12	POS	POS	POS	POS	POS
Animal 2	NS	NS	NS	NS	NS	NS
Animal 3	NS	NS	NS	NS	NS	NS
Animal 4	NS	NS	NS	NS	NS	NS

Evaluation of LFAs on animals infected with Kenya1033 (genotype IX) ASFV strain

Table 4: Summary results from animals infected with Kenya1033 (genotype IX) ASFV strain.

Virus quantity (HAD ₅₀)	Clinical scoring	Assay A	Assay B	Assay C	Assay D	Assay E
Day 0						
Animal 1	0	NEG	NEG	NEG	NEG	NEG
Animal 2	1	NEG	NEG	NEG	NEG	NEG
Animal 3	1	NEG	NEG	NEG	NEG	NEG
Animal 4	1	NEG	NEG	NEG	NEG	NEG
Day 4						
Animal 1	9	NEG	NEG	NEG	NEG	NEG
Animal 2	5	NEG	NEG	NEG	POS	NEG
Animal 3	7	POS	POS	NEG	NEG	NEG
Animal 4	14	NEG	POS	POS	NEG	NEG
Day 5 (PM)						
Animal 1	14	POS	POS	NEG	NEG	NEG
Animal 2	8	NS	NS	NS	NS	NS
Animal 3	14	POS	POS	POS	NEG	NEG
Animal 4	14	POS	POS	POS	NEG	POS
Day 8 (PM)						
Animal 1	NS	NS	NS	NS	NS	NS
Animal 2	17	POS	POS	NEG	NEG	POS
Animal 3	NS	NS	NS	NS	NS	NS
Animal 4	NS	NS	NS	NS	NS	NS

Evaluation of LFAs on animals infected with Uganda1965 (genotype X) ASFV strain

Table 5: Summary results from animals infected with Ken05/Tk1 (genotype X) ASFV strain.

Virus quantity (HAD ₅₀)	Clinical scoring	Assay A	Assay B	Assay C	Assay D	Assay E
Day 0						
Animal 1	1	NEG	NEG	NEG	NEG	NEG
Animal 2	2	NEG	NEG	NEG	NEG	NEG
Animal 3	1	NEG	NEG	NEG	NEG	NEG
Animal 4	1	NEG	NEG	NEG	NEG	NEG
Day 4						
Animal 1	8	NEG	POS	NEG	NEG	NEG
Animal 2	13	POS	POS	NEG	NEG	NEG
Animal 3	13	POS	POS	POS	POS	NEG
Animal 4	12	POS	POS	NEG	POS	NEG
Day 5 (PM)						
Animal 1	18	POS	POS	POS	NEG	NEG
Animal 2	15	POS	POS	POS	NEG	POS
Animal 3	15	POS	POS	NEG	POS	POS
Animal 4	16	POS	POS	NEG	NEG	NEG

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Table 1. ASFV-Ag-negative samples analysed collected from Spanish farms

INgezim® ASFV CROM Ag	Prototype Ag-LFA		
		Pos	Neg
	Pos	0	3
	Neg	0	157

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Table 2. Other ASFV-Ag-negative samples

INgezim® ASFV CROM Ag	Prototype Ag-LFA		
		Pos	Neg
	Pos	2	10
	Neg	0	138

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



*Paired **blood-serum** samples were also evaluated, confirming the maintained sensitivity regardless the matrix used.*

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Table 3. ASFV-Ab-negative samples analysed


Description of samples	Tested (nº)
Field serum samples (ASFV-free region)	162
TB's Ab positive sera	95
PRRSV's Ab positive sera	26
CSFV's Ab positive sera	15

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
Combo CROM

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Combo Ag/Ab cassette

	ANTIGEN POSITIVE	ANTIGEN NEGATIVE
ANTIBODY NEGATIVE		
ANTIBODY POSITIVE		

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Technical Specialist _ BIOLABTEK CO., LTD